



SP Diagnostics - www.spsupport.com




Online Applications Viewer

Guide to online applications viewer

This guide will demonstrate how to access the online applications viewer within the SPi Suite PC software. The online viewer displays data live from SP Diagnostics' central system so it is always up to date.

Version 1.001


1. Open the SP Diagnostics SPi Suite
2. Select "Applications Viewer" from the main menu



3. The SPi Suite will check that you have an internet connection.



4. The Applications viewer will be displayed



5. Use the drop-down lists to filter the list for a specific Manufacturer and/or Model



PO BOX 3467, Chippenham, Wiltshire, SN15 9FH, UK

Sales & Customer Services: + 44 (0)1225 744 777

Facsimile: + 44 (0)8700 941753

Email: sales@spdiagnostics.com

Email: support@spdiagnostics.com

www.spdiagnostics.com

www.spsupport.com